

**TMZF[®] Alloy:
A Titanium Alloy Optimized
for Orthopaedic Implants**




TRANSLATING ADVANCED METALLURGY TECHNOLOGY...

TMZF[®] alloy is an advanced titanium alloy developed in 1986 through Howmedica Osteonics' innovative research in orthopaedic metallurgy. With a unique composition of titanium, molybdenum, zirconium, and iron, it achieves a superior combination of flexibility, strength, and notch resistance when compared to other alloys used in orthopaedic implants, including Ti-6Al-4V. In laboratory testing with Ti-6Al-4V, TMZF[®] alloy further demonstrates improved wear resistance, reducing the potential for generation of particulate metallic wear debris.¹ ***TMZF[®] alloy is specifically tailored for high performance in orthopaedic applications, optimizing the material properties that are key elements in the comfort of your patients and the long-term clinical success of the implant.***


Increased Flexibility

The modulus of implant materials is a core factor in adequate transfer of stress to the surrounding bone. Decreasing material modulus, thereby increasing elasticity, enhances implant-to-bone stress loading and can minimize bone atrophy due to stress shielding. *TMZF[®] alloy, when compared to other alloys used in orthopaedics, has a modulus closer to that of bone, improving stress transfer and potentially decreasing thigh pain in patients which results from implant tip contact with the cortices.*




Greater Strength

Material elasticity is only useful when combined with suitable material strength. *TMZF[®] alloy offers increased implant flexibility without sacrificing component strength. In fact, TMZF[®] alloy has up to a 20% higher yield strength than Ti-6Al-4V alloy.¹*



Improved Notch Resistance

Improved notch resistance provides better tolerance to surface stress concentrations caused by the rigorous biomechanical forces experienced by an implant over time. *Results for notched samples show that TMZF[®] alloy is far superior to Ti-6Al-4V alloy, with a 47% higher fatigue strength.¹ This notch resistance, coupled with improved wear and abrasion resistance, builds a better foundation for implant longevity.*



With its demonstrated advantages in material properties, TMZF[®] alloy, combined with Howmedica Osteonics' clinically successful implant geometries and coating technologies, takes orthopaedic design to new standards of performance.

... INTO ORTHOPAEDIC SOLUTIONS

Howmedica Osteonics offers many implant choices in TMZF® alloy, with geometries and coating technologies designed to accommodate the specific needs of your patients.

Fit & Fill



Meridian® TMZF® Femoral Component

- Two distal diameters available for each proximal body support a broader patient population
- Roughened titanium plasma spray coating mechanically designed to lock the bone and implant to enhance bone ongrowth and initial stability
- Distal split and grooves are designed to improve patient comfort by further reducing distal stiffness without affecting component strength

Citation® TMZF® Femoral Component

- Enlarged anterior geometry and built-in proximal anteversion offers a highly stable fit in the femoral canal
- Howmedica Osteonics' clinically successful PureFix™ HA coating² on a titanium plasma spray surface provides an optimal surface for bone ongrowth
- Small polished neck offers greater ROM with the potential for reduced polyethylene wear

Anatomic



Tapered



Accolade® TMZF® Femoral Component

- Tapered body geometry designed to create a bone-conserving, medial-lateral taper wedge fit which has been successful in 13-year clinical results^{3,4,5,6}
- Incremental neck lengths and standard/extended offset options provide better intra-operative flexibility for adjusting tissue-tensioning and leg length
- Small, trapezoidal neck results in greater ROM in flexion and internal rotation
- PureFix™ HA coating on a titanium plasma spray surface creates further potential for bone ongrowth

References

- ¹ Wang K, Gustavson L, Dumbleton J. The Characterization of Ti-12Mo-6Zr-2Re — A New Biocompatible Titanium Alloy Developed for Surgical Implants. In: Beta Titanium Alloys in the 1990s, Edited by Eylon D. Boyer RR, Koss DA. **The Minerals, Metals & Society, 1993.**
- ² D'Antonio JA, Capello WN, Manley MT, Geesink R. Hydroxyapatite Femoral Stems for Total Hip Arthroplasty: 10- to 13-Year Follow Up. **Clinical Orthopaedics and Related Research.** 393:101-111, 2001.
- ³ Hozack W, Rothman R, Eng K, Mesa J. Primary Cementless Hip Arthroplasty with a Titanium Plasma Sprayed Prosthesis. **Clinical Orthopaedics and Related Research.** 333: 217-225, 1996.
- ⁴ McLaughlin J, Lee K. Total Hip Arthroplasty in Young Patients: 8-13 Year Results Using an Uncemented Stem. **Clinical Orthopaedics and Related Research.** 373:153-163, 2000.
- ⁵ Sakalkale D, Eng K, Hozack W, Rothman R. Minimum 10-Year Results of a Tapered Cementless Hip Replacement. **Clinical Orthopaedics and Related Research.** 362:138-144, 1999.
- ⁶ Keisu KS, Hozack WJ, et al. Primary Cementless Total Hip Arthroplasty in Octogenarians. **The Journal of Bone and Joint Surgery.** 83A:359-363, 2001.

Also available in TMZF® alloy from Stryker Spine:
Reflex™ Anterior Cervical Spinal Plate



Stryker® is a registered trademark of Stryker Corporation.
TMZF®, Citation® and Meridian® are registered trademarks of Howmedica Osteonics Corp.
Accolade™ and PureFix™ are trademarks of Howmedica Osteonics Corp.
©Howmedica Osteonics 2002 Lit. No. LTMZFB 2M 04/02 MC/OI Printed in U.S.A.

stryker®
Howmedica
OSTEONICS

59 Route 17 South
Allendale, NJ 07401
howost.com